



**IMPORTER/DISTRIBUTOR:**

UNIBUY GMBH SP. Z O.O.  
UL. CHMIELNA 2/31  
00-020 WARSZAWA  
EMAIL: SKLEP@INSTRUMEDIC.PL  
TEL. +48 73 74 88 443

# Video Laryngoscope

### Product Introduction

AUG Video Laryngoscope adopts state-of-art technique, possesses the advantage of 2.0M Pixels full-view camera, High-resolution monitor, Unique anti fog capability (no need to wait for preheating), along with portable ergonomic design, enabling clinicians to improve their chances of first-attempt success while minimizing the injury chances may occur during the intubation process by getting a better view of the glottis structure.

### Features:

- 2.0M Pixels camera, full-view high resolution monitor.
- Ergonomic durable design.
- Easy to Operate.
- Utilize 316 medical stainless steel material and blades can be sterilized or through anti-microbial soaking.
- Affordable prices

### System configuration

3.5" Portable color display  
Various specifications of blades  
DISPOSABLE VIDEO LARYNGOSCOPE BLADE (MACINTOSH 2, 3, 4)  
REUSABLE VIDEO LARYNGOSCOPE BLADE (MACINTOSH 2, 3, 4)  
Size 1 and 5 are optional.  
Rechargeable lithium-ion battery  
Instruction manual  
Quick installation guide  
Warranty card

### Operating the VLR-100 & VLD-150 Display

1. Power Key: Switch on/off the device ;
2. M key: menu key / video play Pause key / Enter / time selection.
3. ↑ key: last image or video / increase value .
4. ↓ key: next image or video / Decrease value .
5. Snap key : short press camera / long press video / return .
6. Camera : for being filmed .

Ref: VLD-150





# AUG MEDICAL Video Laryngoscope

## Application:

Difficult airway intubation (CL from I ~ IV grade), Traumatic airway intubation,  
Intubation for cervical activity  
Limited patients  
Replace routine laryngoscope.  
Use for anesthesia department, emergency department, ICU,  
Cardiology, etc.  
Wild rescue, first aid center, ambulance

## Clinical Advantages:

Avoid the injury intubation may cause to the laryngeal structure by  
getting a optimal view of the laryngeal structure  
Unique anti-fog capacity, no need to preheat, save  
time in emergency intubation operation  
Complete specifications of blades (Pediatric Size/Adult  
Size/Difficult Airway Intubation Size optional)  
Ultra-big capacity of the battery, can be continually used for 200  
minutes.



Ref: VLR-100

